

Gene Name

TP53I11, tumor protein p53 inducible protein 11

Synonyms

PIG11

Reactivity

Human

Host

Rabbit

Application

WB (1:500 - 1:1000)
IF/IC (1:100-1:500)

Clonality

Polyclonal

SwissProt

O14683

Gene ID

9537

Concentration

1 mg/ml

Form

Liquid

Purification

The antibody was purified by immunogen affinity chromatography.

Buffer


Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20°C for 12 months. Avoid repeated freeze / thaw cycles.

LifeSab™ PIG11 Rabbit pAb

Cat. #: AD0111 Size: 50 µl/100 µl



Western blot analysis of PIG11 expression in HuvEc (A) whole cell lysates.