

Rabbit Anti-Human IgG Fc Biotin

Cat. #: AD2047 Size: 100 µl/500 µl

Immunogen

Human IgG Fc

Clonality

Polyclonal

Congujation

Biotin

Applications

ELISA, WB, IH, IC

Concentration

1 mg/mL

Dilution

E (1/5000 - 1/20000),
WB (1/5000 - 1/20000),
IH (1/500 - 1/1000),
IC (1/500 -1/1000)

Form

Liquid

Specificity

By immunoelectrophoresis and ELISA this antibody reacts specifically with Human IgG Fc. No antibody was detected against non immunoglobulin serum proteins.

Buffer

Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.05% Sodium Azide.

Purification

The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.

Storage

Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. avoid freeze/thaw cycle.